

PLUS MEDICAL CO.,LTD.

佛山市南海康得福医疗用品有限公司

Add: Shatou Jiujiang Town, Nanhai District, Foshan City, Guangdong Province , China (528208)

Tel:+86-757 8691 0199 Fax:+86-757 8691 6230

SPECIFICATION SHEET

DATE: 2nd, April, 2020

Face Mask (BS EN 14683-2019 Type II)

Ref: SFMW01

Compositions:

1st Layer - 25gsm Blue PPSB-100% polypropylene(Outer Surface)

2nd Layer - 22gsm White Meltblown-100% polypropylene

3rd Layer - 25gsm White PPSB-100% polypropylene(Inner Surface)


Nose Stripe - Polypropylene

Elastic Ear-loop - Spandex + Nylon

Performance:

Standard	Characteristics	Test method	Unit	Type		
				I	II	IIR
BS EN 14683-2019: Medical Face Masks- Requirements and Test Methods	Bacterial Filtration Efficiency (BFE)	EN 14683 Annex B	%	≥95	≥98	≥98
	Differential Pressure (Delta P)	EN 14683 Annex C	Pa/cm2	<40	<40	<60
	Splash Resistance Pressure	ISO 22609	kPa	Not Required	Not Required	≥16.0
	Microbial Cleanliness (Bioburden)	EN ISO 11737-1	cfu/g	≤30	≤30	≤30

Sketch:



Storage Instruction:

Do not place together with highly oxidizable materials.

Please keep away from sources of ignition.

Store at warehouse with good sanitary and aeration, dry and waterproof condition.

No sun light or UV radiation.